

No. 449—BARNES CLAMP, for compression arthrodesis of the knee, complete with four  $3/16"$  x  $10"$ , Smo Stainless Steel Steinman pins.



No. 428—LOWE-BRECK CARTILAGE KNIFE, Stainless.



No. 429—DOWNING CARTILAGE KNIFE, Stainless.




No. 430—KNEE SCALPEL, Stainless.



No. 432—HUNDLEY KNEE KNIFE with removable blade, Stainless.




No. 431—SMILLIE CARTILAGE KNIFE: center, right or left cuts, Stainless.



No. 427—MARTIN CARTILAGE SCISSORS, plain or serrated, Stainless.



No. 419—WALTON CARTILAGE CLAMP, straight or curved, length:  $4\frac{1}{2}"$  or  $8"$ , Stainless.



No. 419-B—WALTON CARTILAGE CLAMP, angular, shaft length:  $3\frac{1}{2}"$  or  $5\frac{1}{2}"$ , Stainless.


No. 426—SHERRILL KNEE RETRACTOR, blade depth:  $1\frac{1}{2}"$ ,  $2"$  or  $2\frac{3}{8}"$ , Stainless.



No. 404—BLOUNT KNEE RETRACTOR, Stainless.

No. 437—McKEEVER CARTILAGE KNIFE, full or medium curve, Stainless.

ALL PRODUCTS OF WRIGHT MANUFACTURING COMPANY ARE OF HIGHEST QUALITY AND ARE UNCONDITIONALLY GUARANTEED AGAINST DEFECTIVE WORKMANSHIP AND MATERIALS.



No. 563-B—CHARNLEY CLAMP for lower extremities (as illustrated) four pin type, Stainless,  $4\frac{1}{2}"$  adjustment, complete with four  $3/16"$  x  $9"$  Smo Stainless Steel Steinman pins.

No. 563-A—CHARNLEY CLAMP for lower extremities (not shown) Stainless,  $5\frac{1}{4}"$  adjustment, complete with two  $3/16"$  x  $9"$  Smo Stainless Steel Steinman pins.




**WRIGHT**

MANUFACTURING COMPANY

880-882 ADAMS

MEMPHIS, TENN.



A true artist.....takes joy in his work.....

## OSTIC helps you make strong, neat casts every time

Ostic is the all-around plaster bandage that's so popular in teaching centers. It's easy to use. Applies with a creaminess that's gentle to your hands. Sets quickly, takes hard use.

Nothing can take the place of professional experience when it comes to building a neat, strong cast. But a good part of your results can depend upon the plaster bandage you use. That's why so many specialists in orthopedic teaching centers use Ostic plaster bandages.

They know they can count on Ostic for consistently fine results. Ostic is fast. It's strong enough to take most any kind of wear. And Ostic delivers at least 97% of the original plaster to the cast. You can work with it all day long.



Each Ostic plaster bandage comes in its own airtight package. Plaster stays fresh indefinitely. Eliminates any danger of presetting.

Gentle creaminess of Ostic is never hard on your hands. You get a choice of two setting times: Fast (5 to 8 minutes) and Extra Fast (2 to 4 minutes). Call your Curity representative for the complete story of OSTIC Plaster Bandages.

**Curity®**

THE **KENDALL** COMPANY  
BAUER & BLACK DIVISION